


## ThetaMetrisis APPLICATION NOTE #012





### Determination of substance concentration in solutions through absorption measurements



**Goal:** The real time measurement of the absorption & concentration of liquid samples.

**Means & Methods:** An FR-Basic VIS/NIR configured to operate in the 360-1000nm spectral range is combined with a FR-Film/Cuvette holder for the accurate real time monitoring of the absorption or transmittance of liquids in standard quartz cuvettes. The substance used in the present study was Methylene blue (MB), a heterocyclic aromatic chemical compound, and the absorbance of its aqueous solutions was measured. At room temperature, MB appears as a solid, odorless, dark green powder that yields a blue solution when dissolved in water. Methylene blue aqueous solutions show absorbance in visible region due to involving of  $\pi$ - $\pi^*$  and  $n$ - $\pi^*$  transitions present in molecules. The specifics of absorption depend on a number of factors, including protonation, adsorption to other materials, and metachromacy - the formation of dimers and higher-order aggregates depending on concentration and other interactions. More often the solvatochromic effect is dominant and dependant on the polarity of solvent.

**Results:** The absorption spectra as they recorded through FR-Monitor are illustrated in the figure below for various concentrations. The concentration range is following the Beer-Lambert law where Absorbance is a linear function of the molar concentration ( $A=\epsilon bC$ ). Thus, calculation of the molar absorptivity (or molar absorption coefficient,  $\epsilon$ ) of the examined chemical compound is attainable. The absorption measurements are performed at 668nm according to bibliography and the results are in good agreement with the literature.



Measurement Set-up. Absorbance calibration curve of Methylene Blue at 668nm at 25°C.

#### Conclusions:

Determination of  $A_{\max}$  wavelength.

Determination of concentration via Beer-Lambert law with high sensitivity and Limit of Detection (LOD), i.e. for Methylene Blue aqueous solutions,  $LOD=4.67 \times 10^{-7} M$ .

The described procedure is reliable, very simple and conveniently applicable in most laboratories with low time of analysis.



[1] T. Mahmood et al, European Academic Research 1(6), 2013, p: 1100, [2] J. Cenens et al, *Clay and Clay Minerals* 36 (3), 1988, p: 214. [3] A. Ghanadzadeh Gilani et al, *Journal of Molecular Liquids* 179, 2013, p: 118.