


## ThetaMetrisis APPLICATION NOTE #005


### Thickness measurement of stack of films by White Light Reflectance Spectroscopy (WLRs)




**Goal:** The accurate measurement of the thickness of films in a stack consisting of 4 layers on a substrate.

**Means & Methods:** WLRs is introduced for the measurement of film thicknesses in the case of stack of dielectric and semiconductor films. All measurements were performed with an FR-Basic in the 400-1000nm spectral regime. The reflection probe had an active spot size of 1mm in diameter. The sample was a Si wafer with areas consisting of  $\text{Si}_3\text{N}_4/\text{SiO}_2$ , poly-Si/ $\text{Si}_3\text{N}_4/\text{SiO}_2$ , PMMA/poly-Si/ $\text{Si}_3\text{N}_4/\text{SiO}_2$ . In all measurements a highly reflective Al coated mirror (NT01-913-533, Edmund Optics) was used for the reference measurements.


**Results:** In fig. 1 the reflectance spectrum from a  $\text{Si}_3\text{N}_4/\text{SiO}_2$  stack on Si wafer is shown. The calculated thicknesses were: 576.63nm for  $\text{SiO}_2$  and 127.78nm for the  $\text{Si}_3\text{N}_4$ . In fig. 2 the reflectance spectrum from a poly-Si/ $\text{Si}_3\text{N}_4/\text{SiO}_2$  stack on Si is shown. The calculated thicknesses were: 578.45nm for  $\text{SiO}_2$  and 144.77nm for the  $\text{Si}_3\text{N}_4$  and 114.82nm for the polycrystalline Si. In fig. 3 the reflectance spectrum from a PMMA/poly-Si/ $\text{Si}_3\text{N}_4/\text{SiO}_2$  stack (4 layers) on Si is shown. The calculated thicknesses were: 575.22nm for  $\text{SiO}_2$  and 147.57nm for the  $\text{Si}_3\text{N}_4$ , 115.73nm for the polycrystalline Si and 132.38nm for PMMA. The individual thickness of each film is calculated simultaneously without any previous knowledge of the film thickness in all cases. In all cases, calculated film thicknesses were in very good agreement with the values obtained by spectroscopic ellipsometry (<1% difference).



**Figure 1:** Typical reflectance spectrum from a  $\text{Si}_3\text{N}_4/\text{SiO}_2/\text{Si}$  sample (raw data without normalization). Black line: sample, Red line: fitting



**Figure 2:** Typical reflectance spectrum from a poly-Si/ $\text{Si}_3\text{N}_4/\text{SiO}_2/\text{Si}$  sample (raw data without normalization). Black line: sample, Red line: fitting



**Figure 3:** Typical reflectance spectrum from a PMMA/poly-Si/ $\text{Si}_3\text{N}_4/\text{SiO}_2/\text{Si}$  sample (raw data without normalization). Black line: sample, Red line: fitting

**Conclusions:** The thickness of at least 4 layers can be calculated simultaneously.